MICROORGANISMOS INDICADORES

Milena Gómez

Sales Application Specialist

milena.gomez@merckgroup.com

MERCK

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1

Generalidades

2

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
3

Mc Media PAD

4

Simplite

Generalidades

MERCK

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CERTIFICATE OF COMPLIANCE LLOYD'S REGISTER QUALITY ASSURANCE

hereby declares that the certification assessment has demonstrated that

MC-Media Pad™ ACplus

Manufactured and supplied by:
JNC Corporation
2-2-1 Otemachi
Chiyoda-ku
Tokyo
Japan, 100-8105

Production site:
JNC Corporation
Yokohama Research Center
5-1, Ookawa
Kanazawa-ku
Yokohama
Japan, 236-8605

has been validated and revealed to be at least equivalent to the reference method as demonstrated by the validation study report. The summary of the validation report is available on the MicroVal website: www.microval.org

Reference method: ISO 4833-1:2013 Microbiology of the food chain — Horizontal method for the enumeration of microorganisms Part 1: Colony count at 30 degrees C by the pour plate technique

Scope: A broad range of foods

The validation and certification has been performed in accordance with EN ISO 16140-2: 2016 and the MicroVal Rules and Certification Scheme version 8.

This certificate is valid only in association with the certificate schedule bearing the same number on which the distributor applicable to this approval is listed.

Certificate no.: 2015LR52

First approval date: 29 January 2018
Expiry date: 28 January 2022

CERTIFICATION

AOAC® Performance TestedSM

Certificate No.
070901

The AOAC Research Institute hereby certifies that the performance of the test kit known as:

MC-Media Pad *E. coli* & Coliform Sheet Medium

manufactured by
JNC Corporation
2-1, Otenachi 2-chome
Chiyoda-ku Tokyo 100-8105
Japan

This method has been evaluated in the AOAC® Performance Tested MethodsSM Program, and found to perform as stated by the manufacturer contingent to the comments contained in the manuscript. This certificate means that an AOAC® Certification Mark License Agreement has been executed which authorizes the manufacturer to display the AOAC Performance TestedSM certification mark along with the statement - "THIS METHOD'S PERFORMANCE WAS REVIEWED BY AOAC RESEARCH INSTITUTE AND WAS FOUND TO PERFORM TO THE MANUFACTURER'S SPECIFICATIONS" - on the above mentioned method for a period of one calendar year from the date of this certificate (June 23, 2017 - December 31, 2017). Renewal may be granted at the end of one year under the rules stated in the licensing agreement.

Deborah McKenzie

Deborah McKenzie, Senior Director
Signature for AOAC Research Institute

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


WORKFLOW

Certificate of Analysis

1.32302.0001 MC-Media Pad Rapid Aerobic Count Convenient culture media for rapid enumeration of aerobic microbial contamination in food and beverages
 Batch 20171221RA

Spec. Values	Batch Values
Inoculum on reference medium (Bacillus subtilis ATCC 6633 (WDCM 00003))	50
Inoculum on reference medium (Staphylococcus aureus ATCC 6538 (WDCM 00032))	64
Inoculum on reference medium (Escherichia coli ATCC 8739 (WDCM 00012))	59
Inoculum on reference medium (Pseudomonas aeruginosa ATCC 27853 (WDCM 00025))	13
Colony count (Bacillus subtilis ATCC 6633 (WDCM 00003))	106
Colony count (Staphylococcus aureus ATCC 6538 (WDCM 00032))	59
Colony count (Escherichia coli ATCC 8739 (WDCM 00012))	69
Colony count (Pseudomonas aeruginosa ATCC 27853 (WDCM 00025))	18
Recovery on test medium (Bacillus subtilis ATCC 6633 (WDCM 00003))	≥ 70 % 212
Recovery on test medium (Staphylococcus aureus ATCC 6538 (WDCM 00032))	≥ 70 % 92
Recovery on test medium (Escherichia coli ATCC 8739 (WDCM 00012))	≥ 70 % 117

Colony count (Pseudomonas aeruginosa ATCC 27853 (WDCM 00025))
 Recovery on test medium (Bacillus subtilis ATCC 6633 (WDCM 00003)) ≥ 70 %
 Recovery on test medium (Staphylococcus aureus ATCC 6538 (WDCM 00032)) ≥ 70 %
 Recovery on test medium (Escherichia coli ATCC 8739 (WDCM 00012)) ≥ 70 %
 Recovery on test medium (Pseudomonas aeruginosa ATCC 27853 (WDCM 00025)) ≥ 70 %
 Colonies (Bacillus subtilis ATCC 6633 (WDCM 00003)) red
 Colonies (Staphylococcus aureus ATCC 6538 (WDCM 00032)) red
 Colonies (Escherichia coli ATCC 8739 (WDCM 00012)) red
 Colonies (Pseudomonas aeruginosa ATCC 27853 (WDCM 00025)) yellowish-brown to red

18
 212 %
 92 %
 117 %
 138 %
 passes test
 passes test
 passes test
 passes test

Incubation: 24 hours at 35 °C, aerobic

Strain selection in accordance with the current version of DIN EN ISO 11133. Test method according to independent validation.

Reference media: Tryptic Soy Agar

Date of release (DD.MM.YYYY) 07.03.2018

Expiry date (DD.MM.YYYY) 21.01.2020

Stefanie Fischer
 Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.


Merck KGaA, Frankfurter Straße 250, 64293 Darmstadt (Germany): +49 6151 72-0
 EMD Millipore Corporation - a subsidiary of Merck KGaA, Darmstadt, Germany
 290 Concord Road, Billerica, MA 01821, USA, Phone: (978) 715-4321
 SALSA Version 692173 /990000516743/ Date: 07.03.2018

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