
MERCK

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Medios de cultivo deshidratados	Art.No.	I&D	Etiqueta	Criterio de prueba	CoA	BSE-Cert.	Hand book	Disponibilidad
Agar Base Cetrimida	1.05284	✓	✓	✓	✓	✓	✓	✓
Agar Sangre Columbia	1.10455	✓	✓	✓	✓	✓	✓	✓
Caldo EE MOSSEL (Harm.)	1.05403	✓	✓	✓	✓	✓	✓	✓
Agar MacCONKEY	1.05465	✓	✓	✓	✓	✓	✓	✓
Caldo MacCONKEY	1.05396	✓	✓	✓	✓	✓	✓	✓
Agar Manitol Sal	1.05404	✓	✓	✓	✓	✓	✓	✓
RCM	1.05411	✓	✓	✓	✓	✓	✓	✓
RVS (Harm.EP)	1.07666	✓	✓	✓	✓	✓	✓	✓
Caldo SAB. 2%-Glucosa	1.08339	✓	✓	✓	✓	✓	✓	✓
Agar SAB. 4%-Glucosa	1.05438	✓	✓	✓	✓	✓	✓	✓
Agar VRBD	1.10275	✓	✓	✓	✓	✓	✓	✓
Agar XLD (Harm.EP)	1.05290	✓	✓	✓	✓	✓	✓	✓
Agar Patata glucosa	1.10130	✓	✓	✓	✓	✓	✓	✓
Caldo NaCl-Peptona	1.10582	✓	✓	✓	✓	✓	✓	✓
Agar Trypticasa Soya	1.05458	✓	✓	✓	✓	✓	✓	✓
Caldo Trypticasa Soya	1.05459	✓	✓	✓	✓	✓	✓	✓

Page

**USP <60>  
TESTS FOR  
*BURKHOLDERIA CEPACIA*  
COMPLEX (Bcc)**

**Millipore®**

Preparation, Separation,  
Filtration & Monitoring Products

# USP <60> TESTS FOR *BURKHOLDERIA CEPACIA* COMPLEX (Bcc)

## <60> MICROBIOLOGICAL EXAMINATION OF NONSTERILE PRODUCTS — TESTS FOR *BURKHOLDERIA CEPACIA* COMPLEX

Printed on: Fri May 07 2021, 07:43:04 AM (EST) Official Status: Currently Official on 07-May-2021 DocId: 1\_GUID-4C3D84A8-B940-42F0-ABC3-AE9128AA8D95\_2\_en-US  
Printed by: Barbara Gerten Official Date: Official as of 01-Dec-2019 Document Type: GENERAL CHAPTER @2021 USPC  
1

**Add the following:**  
▲<60> MICROBIOLOGICAL EXAMINATION OF NONSTERILE PRODUCTS  
—TESTS FOR *BURKHOLDERIA CEPACIA* COMPLEX

**INTRODUCTION**

The tests described in this chapter will allow determination of the absence of *Burkholderia cepacia* complex (Bcc), which can be detected under the conditions described.  
The tests are designed to determine whether a substance or preparation complies with an established specification for microbiological quality and/or to evaluate whether products—especially those for inhalation use or aqueous preparations for oral, oromucosal, cutaneous, or nasal use—contain members of the Bcc.

**GROWTH-PROMOTING AND INHIBITORY PROPERTIES OF THE MEDIA AND SUITABILITY OF TESTS FOR ABSENCE OF Bcc**

Test each batch of ready-prepared medium and each batch of medium prepared from either dehydrated medium or ingredients.

**Preparation of Test Strains**

Use standardized stable suspensions of test strains (see Table 1) NMT 5 passages removed from the original strain culture.

**Table 1. Test Strains of Microorganisms for Growth Promotion and Suitability Testing**

Microorganism	Standard Strain
<i>Burkholderia cepacia</i>	ATCC 25416, NCTC 10743, or CIP 80.24
<i>Burkholderia cenocepacia</i>	ATCC BAA-245 or LMG 16656
<i>Burkholderia multivorans</i>	ATCC BAA-247, LMG 13010, CCUG 34080, CIP 105495, DSM 13243, or NCTC 13007
<i>Pseudomonas aeruginosa</i>	ATCC 9027, NCIMB 8626, CIP 82.118, or NBRC 13275
<i>Staphylococcus aureus</i>	ATCC 6538, NCIMB 9518, CIP 4.83, or NBRC 13276

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Test each batch of ready-prepared medium and each batch of medium prepared from either dehydrated medium or ingredients.

### Preparation of Test Strains

Use standardized stable suspensions of test strains (see *Table 1*) NMT 5 passages removed from the original strain culture.

**Table 1. Test Strains of Microorganisms for Growth Promotion and Suitability Testing**

Microorganism	Standard Strain
<i>Burkholderia cepacia</i>	ATCC 25416, NCTC 10743, or CIP 80.24
<i>Burkholderia cenocepacia</i>	ATCC BAA-245 or LMG 16656
<i>Burkholderia multivorans</i>	ATCC BAA-247, LMG 13010, CCUG 34080, CIP 105495, DSM 13243, or NCTC 13007
<i>Pseudomonas aeruginosa</i>	ATCC 9027, NCIMB 8626, CIP 82.118, or NBRC 13275
<i>Staphylococcus aureus</i>	ATCC 6538, NCIMB 9518, CIP 4.83, or NBRC 13276

# USP <60> TESTS FOR *BURKHOLDERIA CEPACIA* COMPLEX (Bcc)

## <60> MICROBIOLOGICAL EXAMINATION OF NONSTERILE PRODUCTS – TESTS FOR *BURKHOLDERIA CEPACIA* COMPLEX

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**Table 2. Microorganisms for the Growth-Promoting, Inhibitory, and Indicative Properties of the Media**

Medium	Property	Microorganism
<i>Burkholderia cepacia</i> selective agar	Growth-Promoting and Indicative	<i>Burkholderia cepacia</i> , <i>Burkholderia cenocepacia</i> , or <i>Burkholderia multivorans</i>
	Inhibitory	<i>Pseudomonas aeruginosa</i> , <i>Staphylococcus aureus</i>

# USP <60> TESTS FOR *BURKHOLDERIA CEPACIA* COMPLEX (BCc)

## <60> MICROBIOLOGICAL EXAMINATION OF NONSTERILE PRODUCTS – TESTS FOR *BURKHOLDERIA CEPACIA* COMPLEX

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#### *Burkholderia cepacia* Selective Agar

Prepare BCSA as directed in Table 5. When preparing media in-house, first prepare the base ingredients without the antibiotics. Adjust the pH so that after sterilization it is  $6.8 \pm 0.3$  at 25°. Sterilize in an autoclave using a validated cycle. Cool the base medium to 45°–50° and add a 1% solution of the sterile filtered antibiotics, mix, and pour into the plates.

Table 5

Casein peptone	10.0 g
Lactose	10.0 g
Sucrose	10.0 g
Sodium chloride	5.0 g
Yeast extract	1.5 g
Phenol red	0.08 g
Gentamicin	10.0 mg
Vancomycin	2.5 mg
Crystal violet	2.0 mg
Polymyxin B	600,000 U
Agar	14.0 g
Deminerlized water	1000 mL ▲ (USP 1-Dec-2019)

# USP <60> TESTS FOR *BURKHOLDERIA CEPACIA* COMPLEX (Bcc)

## <60> MICROBIOLOGICAL EXAMINATION OF NONSTERILE PRODUCTS

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**CAS Number:**  
**Formula:**  
**Formula Weight:**  
**Expiration Date:** JAN 2024  
**Quality Release Date:** 16 MAR 2021

TEST	SPECIFICATION	RESULT
APPEARANCE (COLOR)	FAINT RED TO LIGHT RED AND PINK	FAINT RED
APPEARANCE (FORM)	POWDER	POWDER
MICROBIOLOGY TESTS	PASS AS DESCRIBED IN DATA SHEET	PASS
SOLUBILITY (COLOR)	COLORLESS TO DARK ORANGE AND COLORLESS TO DARK RED	DARK ORANGE
SOLUBILITY (TURBIDITY)	CLEAR	CLEAR
PH	6.85 ± 0.2 AT 25 C	6.9
GEL STRENGTH	CORRESPONDS TO REQUIREMENTS	CORRESPONDS
SOLUBILITY (METHOD)	51.6 G IN 1000 ML WATER	51.6 G IN 1000 ML WATER

# Complejo *Burkholderia cepacia*

## CoA

Cultural characteristics after 48-72 hours at 30-35°C. Inoculum: <100 CFU (Productivity) according USP, practical 50-100 CFU according to ISO 11133:2014/Amd1:2018/Amd2:2020. Spiral Plate Method.

Organisms (ATCC®, <b>WDCM</b> )	Recovery	Remarks
<i>Burkholderia cepacia</i> (25608)	≥50%	Greenish-brown colonies with yellow halo
<i>Burkholderia cepacia</i> (25416*)	≥50%	Greenish-brown colonies with yellow halo
<i>Burkholderia cenocepacia</i> (BAA-245*)	≥50%	White colonies
<i>Burkholderia multivorans</i> (BAA-247*)	≥50%	White colonies
<i>Staphylococcus aureus</i> (6538*/ <b>00032</b> )	inhibited	w. selective supplement
<i>Pseudomonas aeruginosa</i> (9027*/ <b>00026</b> )	Inhibited (partial to complete)	w. selective supplement

\*all test strains recommended by USP-NF chapter 60